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
Read these Instructions carefully
before use.

BlueShot II

Dental Light Curing Unit

Instructions for Use

SHOFU INC.

Introduction

Thank you for purchasing our SHOFU BlueShot II. The instructions mentioned below provide necessary information on successful operation, daily maintenance and safe use of this unit.

Please read the instructions carefully before use to attain maximum benefits and ensure durability of this unit. Kindly keep this manual for your future reference.

Notice

- **No part of this manual may be copied or reproduced in any form by any means.**
 - **Some information in this manual may become inconsistent with the product due to product software updates.**
 - **The information in this manual is subject to change without prior notice. All relevant issues have been considered in the preparation of this document. However, should you notice an omission or any questionable item in this document, please feel free to notify SHOFU.**
 - **In case of missing pages or disorder, please contact our local dealer for replacement.**
 - **SHOFU will repair and replace main body and its system in case of trouble within the warranty period. However, we do not assume any responsibility for subsequent damages arising from troubles during operation.**
-

Contents

Introduction.....	ii
Notice	ii
Contents	iii
Indications	iv

1 Safety Information.....	1
Symbols.....	1
Notes and Signs	1
Installation and Connection	1
Operation.....	2
Maintenance	3

2 List of Components	4
-----------------------------------	----------

3 Parts and Functions.....	5
Symbols on Product.....	6

4 Installation and Connection	7
--	----------

5 Directions for Use	9
-----------------------------------	----------

6 Selection of Curing Mode.....	12
--	-----------

7 Durable Years	13
------------------------------	-----------

8 Maintenance	13
Cleaning and Sterilization	13
Maintenance and Replacement of the Battery	13

9 Troubleshooting	14
--------------------------------	-----------

10 Storage and Transport.....	15
--------------------------------------	-----------

11 Disposal	15
--------------------------	-----------

12 Specification.....	16
------------------------------	-----------

13 Accessories	17
-----------------------------	-----------

14 Warranty.....	17
-------------------------	-----------

15 EMC (Electromagnetic compatibility).....	17
--	-----------

Indications

For polymerization of light curing materials including the light cure catalyst initiated by the light with wavelength range 440 nm to 490 nm.

1 Safety Information

It is important to understand the following information to use this product safely.

Symbols

 **Warning:** Indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

 **Caution:** Indicates a potentially hazardous situation which, if not avoided, may result in minor/moderate injury or damage of property.

Notes and Signs

Other notes and signs used in this manual are shown below.

Note • Indicates that potential malfunctions may occur if the description is ignored.

Installation and Connection

Warning

- Keep away from sources of ignition to avoid fire or explosion.
 - Protect the Handpiece, Charger Base and Battery Pack from any liquid to avoid electric shock or functional error.
-

Caution

- Securely connect the dust-free power plug to the outlet. Use this unit in well ventilated areas. Overheating may cause fire.
 - Place the Charger Base on flat and stable surface not to get injured by the dropping unit.
-

Warning

- Do not plug the power cord with wet hands. This may result in an electric shock.
 - In case of fluid leak from battery or partial discoloration/deformation or nasty smell of the Handpiece, immediately discontinues use and send back to your dealer for repair.
 - In case of partial discoloration, deformation or nasty smell of the Charger Base, immediately disconnect the power plug and repair the base.
 - In case battery fluid contacts the eye, immediately flush with copious amount of clean water and seek medical attention.
 - Wear the protest glasses or the Anti-glare shield to avoid direct eye contact with LED light and also protect patient's eyes during exposure. And avoid to see the reflected LED light directly. As it may affect eyesight.
 - Do not use this product on patients with a history of cataract, retinal damage or allergic reaction to light.
 - Do not use this product on a person with a pacemaker as it may cause malfunction of the pacemaker.
 - Avoid exposing the light on intraoral soft tissue as it may affect the soft tissue.
-

Caution

- Disconnect the power cord by grasping the power plug. This may cause injuries or burns.
 - Avoid damage or use of unnecessary force on the power cord. This may cause an electric shock or fire.
 - Do not exposure excessive time on same place as it may cause damage to gingiva and pulp.
 - Should battery fluid contact the skin or clothes, immediately wash it off with copious amount of water to avoid skin damage.
 - Take the battery out of the Handpiece if it won't be used for more than 1 week.
 - Connect the power plug only when battery is being charged. Insulation degradation may cause electric shock or fire.
 - Do not use this product for any purposes other than dental use.
-

Maintenance

Warning

- Use only the battery supplied to avoid fluid leak or explosion.
 - Never disassemble or overhaul the product. Unusual performance may cause injuries, electric shock, or fire.
-

Caution


- If LCD display malfunctions, remove the battery from the Handpiece and reinsert it.
 - Users must be in charge of operation, maintenance and check.
-

2 List of Components


Please check the contents described below when removing the contents from the original packaging.


Handpiece


Charger Base


AC Adapter

Battery Pack

Light Probe
(ϕ 8 mm)


Anti-glare Shield

- Manual (Instructions for Use)
- Notes on use of medical electrical equipment (for Japan)
- Warranty cards
- Japanese Package Insert
- Label

3 Parts and Functions

1 Handpiece

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For selection of Turbo Mode. Each time the button is pressed, the exposure time changes to 5 sec., 10 sec., and 20 sec.

1-B Time Button

For selection of exposure time. Each time the button is pressed, the exposure time changes on each mode as follows.

- Turbo Mode : 5 sec., 10 sec., 20 sec.
- High Mode : 10 sec., 20 sec., 30 sec.
- Soft Start Mode : 20 sec.
- Pulse Mode : 20 sec., 30 sec., CON (continuous)
- Normal Mode : 20 sec., 30 sec., CON (continuous)

* Maximum exposure time of CON is approximately 3 min.

To stop the light exposure, press any button on the Handpiece.

1-C Mode Button

For selection of exposure mode.

High Mode → Soft Start Mode → Pulse Mode → Normal Mode → High Mode

1-D Start Button (Power ON/OFF Button)

To start the light exposure, press the button.

To turn ON, press the button when display is disappeared.

To turn OFF, keep pressing the button when the Handpiece is ON.

* To stop the light exposure, press any button on the Handpiece.

1-E LCD Display

Selected mode, exposure time and remaining battery are displayed on the LCD Display in Graph Mode or Text Mode.

2 Light Probe

Effectively exposures the LED light.

3 Anti-glare Shield

It is placed on the Light Probe to protect eyes from exposure light.

4 Battery Pack

4-A Charging terminal area

4-B Connecting terminal area

5 Charger Base

5-A Light display to indicate the charging process with green and red lights.

5-B Inlet

6 AC Adapter

Symbols on Product

	Trade Mark
BlueShot II	Product Logo Mark
	TYPE BF APPLIED PART
	CLASS II equipment
 Li-ion	Recycle Lithium-ion battery

4 Installation and Connection

- Refer to 1 Safety Information about Installation and Connection.


1. Connection between AC Adapter and Charger Base

- (1) Securely insert the connector of the AC Adaptor into the Inlet at the back of the Charger Base.
- (2) The AC adaptor is compatible with a power outlet of 100 V to 240 V.


Note

- Use only the original AC adaptor supplied.
- When the power plug is connected, both red and green lights flash instantly which is not a malfunction or an alarm sign.


2. Inserting Battery Pack

- (1) Carefully insert the Battery Pack into the Handpiece as directed by the arrow. Avoid forcing the Battery Pack.


3. Inserting Light Probe

- (1) Slowly insert the Light Probe into the Handpiece. Push the Light Probe in till it firmly sets into position.

4. Placing Anti-glare Shield

- (1) Place the Anti-glare Shield on the Light Probe.

5. Charging Handpiece


- (1) Insert the Handpiece into the vertical slot at the front of the Charger Base.
The charger light at the front of the Charger Base will appear.

Red light On : Charging in progress
Green light On : Charging complete

- If the battery power level is extremely low, red light may flicker while recharging, which is not a malfunction.
- It will take approximately 3 hours before the battery will be fully charged

Note Depending on the volume of the remaining battery, the sign of battery indicator varies.


5 Directions for Use

- Refer to 1 Safety Information about Operation.

1. Activating Handpiece

- (1) Press the Start Button, LCD Display indicates standby mode and then the mode used last with a short audible tone.
- (2) Select the curing mode and time by pressing the Turbo Mode Button, Mode Button and Time Button.

Initial Screen


Mode used last

2. Selection of the display mode


Keep pressing Mode Button to select the display mode (Graph or Text).


3. LCD Display

Graph Mode and Text Mode are available.

<Graph Mode>


<Text Mode>


Mode	Abbreviation of curing mode
Turbo	Tur.
High	High
Soft Start	Soft
Pulse	Pul.
Normal	Norm

Graph Mode Display

- (1) Abbreviation of curing mode
- (3) Remaining exposure time
- (4) Battery power level
- (5) Output graph

Text Mode Display

- (1) Abbreviation of curing mode
- (2) Selected exposure time
- (3) Remaining exposure time
- (4) Battery power level

Note

Initially, the selected exposure time is displayed and subsequently the remaining exposure time is displayed during the curing process.

4. Start of exposure

- (1) Press Start Button to start exposure with the selected mode.
- (2) Press any button to stop the light exposure.

5. Sleep Mode

If unit is left unused for more than 5 minutes, all functions are automatically turned off. To turn off the device anytime, keep pressing Start Button.

6. Set-up memory


The unit records the setting used last.

6 Selection of Curing Mode

Curing mode can be selected by pressing the Mode Button or Turbo Mode Button. Five curing programs are available.


1. Turbo Mode

The output light intensity is about 1,400 mW/cm². The exposure time can be 5 sec., 10 sec. and 20 sec.


2. High Mode

The output light intensity is about 1,000 mW/cm². The exposure time can be 10 sec., 20 sec. and 30 sec.


3. Soft Start Mode

In the 5 seconds, output light intensity is increasing from 0 mW/cm² to about 650 mW/cm² and then increases to about 1,000 mW/cm² during the next 15 seconds. The exposure time can be 20 sec.


4. Pulse Mode

The output light intensity starts at about 1,000 mW/cm² and holds for 0.1 seconds followed by a rest period of 0.05 seconds in repetition. The exposure time can be 20 sec., 30 sec. and CON (max. 3 min.).

5. Normal Mode

Its output light intensity is about 650 mW/cm². The exposure time can be 20 sec., 30 sec. and CON (max. 3 min.).

7 Durable Years

Durable year are 7 years from the date of purchase under the unit is kept by 8 Maintenance.

8 Maintenance

- Refer to 1 Safety Information about Maintenance.

Cleaning and Sterilization

- Avoid wetting the Handpiece and Charger Base to prevent electrical shock. Ensure that the Charger Base and Battery Pack are always separated from the Handpiece prior to cleaning.
- Clean the Handpiece and Charger Base gently with a clean, soft, lint-free cloth. Avoid using solvents containing oil or more than 50% of ethanol or acetone, sharp or pointed tools to clean as they may damage or scratch the plastic.
- Clean the Light Probe with a soft cloth and disinfect. The Light Probe is autoclavable less than 135°C.
- Remove restorative materials stuck to the Light Probe gently with a plastic spatula. Avoid using sharp tools as they may scratch the surface of the Light Probe and affect the intensity of the light.
- Avoid using alcohol and solvents to clean Anti-glare Shield.

Maintenance and Replacement of the Battery

- Exposure to high room temperature may shorten the battery life. The recommended room temperature ranges between 15 to 25°C. Extremes of room temperature will affect the operation of the LED light.
- It is recommended for the rechargeable battery to be charged immediately before use to maximize the efficiency.
- Remove the battery from the Handpiece, if not being used for more than one week.
- The rechargeable battery operating under normal conditions can be used for around 400 charging cycles (This may vary depending on the using condition).
- If the running time declines with a full-charged battery, replace the battery before use.

9 Troubleshooting

If your curing light unit is not operating normally, discontinue use and send back to your dealer for inspection and repair.

If your unit corresponds following Troubles, the unit may be normal. Check your unit again along Cause and Measure before return the unit.

No.	Trouble	Cause and Measure
1	Control button is not operating normally.	<ul style="list-style-type: none">• If the control button is pressed with excess force, the power will automatically be switched off. In such a case press the button gently.• In case it still does not operate normally, reset the equipment and operate after charging.
2	Display screen does not light.	<ul style="list-style-type: none">• If the battery is completely discharged or connected incorrectly, display will not appear.• Recharge the battery or reconnect the battery to reset the device.
3	Power off while in operation.	<ul style="list-style-type: none">• If the Start Button is pressed for more than 2 sec., power will be switched off. Press the Start Button again to restart the device.• Check the battery connection.
4	Not able to charge.	<ul style="list-style-type: none">• If the power cord is not correctly connected to the Charger Base, or the Handpiece is not properly placed on the Charger Base, the battery can not be recharged.• Clean the Charger Base or the charging terminal.
5	When light intensity is poor.	<ul style="list-style-type: none">• If the Light Probe has dirt or is damaged and the Battery Pack is not charged enough, intensity of light becomes poor.• Replace with new Light Probe or new Battery Pack.

10 Storage and Transport

- Keep this packing materials until disposal to use for transport and storage.
- Send back to your dealer if the unit need repair.
- Only dental professionals must be in charge of storage and management of this unit.
- Keep stable, away from tilt, vibration and impact (during transportation as well).
- Keep away from direct sunlight, contamination and moisture.

11 Disposal

- Do not dispose pollutant including the unit, accessories, packing materials etc.
- Do not dispose the unit and accessories as general waste. Follow each local law and regulation.

12 Specification

Handpiece	Model	BS2-01
	Operating voltage	DC 3.7 V
	Dimensions	L217 × W97 × H38 mm
	Weight	115 g
	Usage environment	Temperature 10°C to 40°C, Humidity 30% to 85%
	Storage environment	Temperature 0°C to 50°C, Humidity 10% to 85%
	Protect against electric shock	TYPE BF APPLIED PART
	Output light intensity	1400 mW/cm ² (Turbo Mode)
	Wavelength range	455 nm to 475 nm
Charger Base	Model	BS2-02
	Operating voltage	DC 5 V
	Dimensions	D179 × W176 × H138 mm
	Weight	240 g
Battery Pack	Model	BS-BP
	Type	Lithium-ion battery
	Nominal voltage	DC 3.7 V
	Rating capacity	2600 mAh
	Dimensions	L108 × H34 mm
	Weight	65 g
Light Probe	Diameter	φ 8 mm
	Length	95 mm
Anti-glare Shield	Material	Acrylic
	Dimensions	L90 × H60 × T3 mm
AC Adaptor	Power supply	AC 100 V to 240 V
	Power output	DC 5 V

* This specification is subject to change without prior notice.

13 Accessories

- Light Probe (PN5112)
- Anti-glare Shield (PN5113)
- Battery Pack (PN5114)
- AC Adapter (PN5115)

14 Warranty

This product is shipped after strict inspection. Should the product have any problem under normal usage within one year from the date of purchase, repair will be free of charge. However, repair of consumables (Light Probe, Battery Pack and Anti-glare Shield etc) are not covered by the warranty. For more information, refer to the Warranty card enclosed.

15 EMC (Electromagnetic compatibility)

Note

- This unit may be affected by portable and mobile type RF telecommunication equipment.
- This unit needs special notice about EMC. Please obey EMC information below.

Guidance and manufacturer's declaration – electromagnetic emissions		
This Model is intended for use in the electromagnetic environment specified below. The customer or the user of this Model should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This Model is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The [EQUIPMENT or SYSTEM] is intended for use in the electromagnetic environment specified below. The customer or the user of the [EQUIPMENT or SYSTEM] should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T (> 95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 sec	< 5% U_T (> 95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the [EQUIPMENT or SYSTEM] requires continued operation during power mains interruptions, it is recommended that the [EQUIPMENT or SYSTEM] be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

This Model is intended for use in the electromagnetic environment specified below. The customer or the user of this Model should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of this Model, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ <p>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this Model is used exceeds the applicable RF compliance level above, the Model should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this Model.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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SHOFU DENTAL ASIA-PACIFIC PTE. LTD. 10 Science Park Road, #03-12, The Alpha, Science Park II, Singapore 117684