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
Read these Instructions carefully
for safe use.

BlueShot

Dental Light Curing Unit

Instructions for Use

SHOFU INC.

Introduction

Thank you for purchasing the Shofu BlueShot. The instructions described below have been designed to provide necessary information on successful operation, daily maintenance and safety for use of this unit.

Please read the Instructions carefully before use to ensure maximum benefits during function and durability of this unit.

Contents

Introduction	ii
Contents	ii
Features	iii
Application	iii
1 Safety Information	1
Symbols	1
Notes and Signs	1
Installation and Connection	2
Operation	3
Maintenance	4
2 List of Components	5
3 Parts and Functions	6
4 Installation and Connection	8
5 Directions for Use	10
6 Selection of Curing Modes	13
7 Maintenance	14
Cleaning and sterilization	14
Maintenance and replacement of the battery	14
8 Trouble Shooting	15
9 Specification	16
10 Accessories	16
11 Warranty	16

Features

- Easy-to-read full color, liquid crystal display.
- 4 types of curing modes.
- Approximately 100 ten-second cures (at High Mode) with fully charged battery.
- Lithium-ion battery that can be recharged any time without depleting the battery.
- Automatic power-off if left unused for more than 4 minutes.
- Autoclavable Light Probe at 135°C (except for the Handpiece) Ergonomic and user-friendly design.

Application

Polymerization of light-cure catalyst containing dental materials with a wavelength spectrum of 440 - 490nm.

1 Safety Information

It is important to understand the following information to use this product safely.

Symbols

 **Warning:** Indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

 **Caution:** Indicates a potentially hazardous situation which, if not avoided, may result in minor / moderate injury or in damage of property.

Notes and Signs

Other notes and signs used in this manual are shown below.

Note • Indicates that potential malfunctions may occur if the description is ignored.

Reference • Supplementary information to facilitate better operation.

Installation and Connection

Warning:

- The adaptor should be connected or disconnected from the Charger Base before connecting the power cable to prevent an electric shock.
 - Keep away from sources of ignition to avoid fire or explosion.
 - Protect the Handpiece, Charger Base and Battery Pack from any liquid to avoid electric shock or functional error.
 - Place the Charger Base on a flat and stable surface. Do not drop the Handpiece and Charger Base to avoid internal defect that may cause overheating or fire.
 - Do not place the electrical cord near a gas burner or sources of heat. Should the cord catch on fire, replace with the new one to avoid short-circuit.
-

Caution:

- Do not use or leave the unit in areas of high temperature or exposed to direct sunlight. Internal circuit defect may cause overheat or fire.
 - Securely connect the dust-free power plug to the outlet. Use this device in well-ventilated areas. Overheating may cause fire.
 - Do not touch the terminal areas of Charger Base or Battery with fingers. Do not insert conductive objects, such as metal wire or safety-pin. It may cause electric shock or fire.
 - Do not use any charger other than supplied. It may cause fire due to overheating.
-

Operation

Warning:

- Do not plug the power cord with wet hands. This may result in an electric shock.
 - In case of fluid leak from battery or partial discoloration / deformation of the Handpiece, immediately discontinue use.
 - In case of partial discoloration or deformation of the Charger Base, immediately disconnect the power plug and repair the base.
 - In case battery fluid contacts the eye, immediately flush with copious amount of clean water and seek medical attention.
 - Avoid direct eye contact with the LED light nor direct the light to someone's eyes as it may affect the eyesight.
 - Do not use this product on patients with a history of cataract, retinal damage or allergic reaction to light.
 - Do not use this product on a person with a pacemaker as it may cause malfunction of the pacemaker.
 - Avoid exposing the light on intraoral soft tissue as it may affect the soft tissue.
-

Caution:

- Disconnect the power cord by grasping the power plug. This may cause injuries or burns.
 - Avoid damage or use of unnecessary force on the power cord. This may cause an electric shock or fire.
 - Should battery fluid contact the skin or clothes, immediately wash it off with copious amount of water to avoid skin damage.
 - Take the battery out of the Handpiece if it won't be used for more than 1 week.
 - Connect the power plug only when battery is being charged. Insulation degradation may cause electric shock or fire.
 - Do not use this product for any purposes other than dental use.
-

Maintenance

Warning:

- Use only the battery supplied to avoid fluid leak or explosion.
 - Never disassemble or overhaul the product. Unusual performance may cause injuries, electric shock, or fire.
-

Caution:


- If the liquid crystal display malfunctions, remove the battery from the Handpiece and reinsert it.
 - Users must be in charge of operation, maintenance and check.
-

2 List of Components


Please check the contents described below when removing the contents from the original packaging.


Handpiece


Charger Base


AC Adapter

Battery Pack

Light Probe
(Ø 8 mm)


Anti-glare Shield

- Manual (Instructions for Use)
- Notes on use of electro-medical apparatus (for Japan)
- Warranty cards
- Japanese DFU
- Label

3 Parts and Functions

1 Handpiece

1-A High Mode Button

For selection of high mode. Each time the button is pressed, the curing time changes to 10 sec., 20 sec., and 30 sec.

1-B Time Button

For selection of exposure time. Each time the button is pressed, the curing time depending on the curing mode.

High Mode : 10 sec., 20 sec., 30 sec.

Normal Mode : 20 sec., 30 sec., CON (continuous)

Soft Start Mode : 20sec.

Pulse Mode : 20 sec., 30 sec., CON (continuous)

* Maximum exposure time of CON is approximately 3 min.

1-C Mode Button

For selection of curing mode.

Normal Mode → Soft Start Mode → Pulse Mode → Normal Mode

1-D Start Button (Power ON/OFF Button)

For switching the power ON/OFF or to start the light exposure.

* To stop the light exposure, press any button on the Handpiece.

1-E Color LCD Display

Selected mode, exposure time and remaining battery are displayed on the LCD display in Graph Mode or Text Mode.

2 Light Probe

Effectively exposures the LED light.

3 Anti-glare Shield

It is placed on the Light Probe to protect eyes from exposure to light.

4 Battery Pack

4-A Charging terminal area

4-B Connecting terminal area

5 Charger Base

5-A Light display to indicate the charging process with green and red lights.

5-B Inlet

6 AC Adapter

4 Installation and Connection

Refer to the information supplied on safety precautions during Installation and Connection.


1. Connection between AC Adapter and Charger Base

- (1) Securely insert the connector of the AC Adaptor into the Inlet at the back of the Charger Base.
- (2) The AC adaptor is compatible with a power outlet of 100V-240V.


Note

- Use only the original AC adaptor supplied.
- When the power plug is connected, both red and green lights flash instantly which is not a malfunction or an alarm sign.


2. Inserting Battery Pack

- (1) Carefully insert the Battery Pack into the Handpiece as directed by the arrow. Avoid forcing the battery pack.


3. Inserting Light Probe

- (1) Slowly insert the Light Probe into the Handpiece. Push the Light Probe in till it firmly sets into position.

4. Placing Anti-glare Shield

- (1) Place the Anti-glare Shield on the Light Probe.

5. Charging Handpiece

- (1) Insert the Handpiece into the vertical slot at the front of the Charger Base. The charger light at the front of the Charger Base will appear.

Red light On: Charging in progress


Green light On: Charging process is complete

- If the remaining battery is extremely low, red light may flickers while recharging, which is not a malfunction.
- It will take approximately 3 hours before the battery will be fully charged

Note


Depending on the volume of the remaining battery, the sign of battery indicator varies.

5 Directions for Use

1. Activating Handpiece

- (1) Press the Start Button gently for 2 sec. After 1 sec. the display screen will be on standby mode and display the mode used last with a short audible tone. (When battery pack is inserted, the unit is automatically switched on.)
- (2) Select the curing mode and time by pressing the High Mode Button, Mode Button and Time Button.

Initial Screen


Mode used last

2. Selection of the display mode


Select the display mode (Graphic or Text) by pressing Mode Button for 1.5 sec. Two short beeps will indicate that the selection is complete.


3. LCD display

Graph Mode and Text Mode are available.

<Graph Mode>

Graph Mode Display


- (1) Abbreviation of curing mode
- (2) Remaining Curing time
- (3) Battery
- (4) Output graph

<Text Mode>

Text Mode Display

- (1) Abbreviation of curing mode
- (2) Remaining curing time
- (3) Battery
- (4) Program Type

- (1) Curing time

Note Initially, the selected curing time is displayed and subsequently the remaining time is displayed during the curing process.

4. Start of exposure


(1) Press Start Button to start exposure with the selected mode.

Press any button if you wish to end the curing process.

5. Battery icon

This icon indicates the power remaining in the battery pack.

Recharge or replace the battery pack when “LOW BATTERY” sign appears on the display screen.

6. Sleep Mode

If unit is left unused for more than 4 minutes, all functions are automatically turned off with a short beep.

Press Start Button for 2 sec. to turn off the device any time.

7. Audible sound

There are 3 types of audible beeps for easy understanding.

- a. Long beep : When putting battery into the handpiece or at the end of light exposure
- b. Two short beeps : When unit is switched on. At the start of light exposure.
- c. Short beep : When unit is swited off , When changing curing mode, Every 10 sec. during exposure (every 5 sec. in High Mode).

8. Set-up memory

The handpiece will record the setting used last.

6 Selection of Curing Mode


Curing mode can be set by pressing the Mode Button or High Mode Button. Four curing programs are available in this LED light.

1. High Mode

This mode provides strong light intensity in a short time.

The output power is about 1000mW/cm².


The exposure time can be 10 sec., 20 sec. and 30 sec.

2. Normal Mode

Its output power is about 650mW/cm².


The exposure time can be 20 sec., 30 sec. and CON (max. 3 min.).

3. Soft Start Mode

In the 5 seconds, output power is raised from 0mW/cm² to 650mW/cm² and then increases to 1000mW/cm² during the next 15 seconds.


The exposure time can be 20 sec.

4. Pulse Mode

The output power starts at 1000mW/cm² and holds for 0.1 seconds followed by a rest period of 0.05 seconds in repetition.

The exposure time can be 20 sec., 30 sec. and CON (max. 3 min.).

7 Maintenance

Cleaning and sterilization

- Avoid wetting the Handpiece and Charger Base to prevent electrical shock. Ensure that the Charger Base and Battery Pack are always separated from the Handpiece prior to cleaning.
- Clean the Handpiece and Charger Base gently with a clean, soft, lint-free cloth. Avoid using solvents containing oil or more than 50% of ethanol or acetone, sharp or pointed tools to clean as they may damage or scratch the plastic.
- Clean the Light Probe with a soft cloth and disinfect. The Light Probe is autoclavable.
- Remove restorative materials stuck to the Light Probe gently with a plastic spatula. Avoid using sharp tools as they may scratch the surface of the Light Probe and affect the intensity of the light.

Maintenance and replacement of the battery

- Exposure to high room temperature may shorten the battery life. The recommended room temperature ranges between 15 to 25°C normal temperature. Extremes of room temperature will affect the operation of the LED light.
- It is recommended for the rechargeable battery to be charged immediately before use to maximize the efficiency.
- Remove the battery from the Handpiece, if not being used for more than one week.
- The rechargeable battery operating under normal conditions can be used for around 400 charging cycles (This may vary depending on the using condition).
- If the running time declines with a full-charged battery, replace the battery before use.
- Exposure to high room temperature may shorten the battery life. The recommended room temperature ranges between 15 to 25°C normal temperature.

8 Trouble Shooting

If your curing light unit is not operating normally, please refer to the following points.

No.	Trouble	Cause and Solution
1	Control button is not operating normally.	<ul style="list-style-type: none">• If the control button is pressed with excess force, the power will automatically be switched off. In such a case press the button gently.• In case it still does not operate normally, reset the equipment and operate after charging.
2	Display screen does not light.	<ul style="list-style-type: none">• If the battery is completely discharged or connected incorrectly, display will not appear.• Recharge the battery or reconnect the battery to reset the device.
3	Power off while in operation.	<ul style="list-style-type: none">• If the Start Button is pressed for more than 2 sec., power will be switched off. Press the Start Button again to restart the device.• Check the battery connection.
4	Not able to charge.	<ul style="list-style-type: none">• If the power cord is not correctly connected to the Charger Base, or the Handpiece is not properly placed on the Charger Base, the battery can not be recharged.• Clean the Charger Base or the charging terminal.
5	When intensity of light is poor.	<ul style="list-style-type: none">• If the Light Probe has dirt or is damaged and the Battery Pack is not charged enough, intensity of light becomes poor. Replace with new Light Probe or new Battery Pack.

9 Specification

Handpiece	Model	BS-01
	Operating voltage	DC3.7V
	Dimensions	Approx. L215 X W96 X H37mm
	Weight	Approx. 175g (including the battery pack)
	Light output	Approx. 1000mW/cm ² (High mode)
	Wavelength range	440-490nm
Charger Base	Model	BS-02
	Operating voltage	DC5V
	Dimensions	Approx. D175 X W178 X H137mm
	Weight	Approx. 225g
Battery Pack	Model	BS-03
	Type	Lithium-ion battery
	Nominal voltage	DC3.7V
	Rating capacity	2400mAh
	Dimensions	Approx. L107 X Max. H33mm
	Weight	Approx.65g
Light Probe	Curing area	Ø 8mm
	Length	95mm
Anti-glare Shield	Material	Acrylic
	Dimensions	Approx. L90 X H60 X T3mm
AC Adaptor	Power supply	AC 100V-240V
	Power output	DC5V

* This specification is subject to change without prior notice.

10 Accessories

- Light Probe (PN5112)
- Anti-glare Shield (PN5113)
- Battery Pack (PN5114)
- AC Adapter (PN5115)

11 Warranty

This product is shipped after strict inspection. Should the product have any problem under normal usage within one year from the date of purchase, repair will be free of charge. However, repair of consumables (Light Probe, Battery Pack and Anti-glare Shield etc) are not covered by the warranty. For more information, refer to the Warranty card enclosed.

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SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan

Distributed by

SHOFU INC., Singapore Branch Blk 215 Henderson Road #03-01, Singapore 159554

SHOFU INC., Shanghai Rep. Office No. 900 Yishan Road, Room 1004-1006, Hi-Tech Building, Shanghai 200233, China

SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA

SHOFU DENTAL PRODUCTS LTD. Duke's Factory, Chiddingstone Causeway, Tonbridge, Kent, TN11 8 JU, UK

EU REPRESENTATIVE: SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany