



VINTAGE AL&ZR



Opal MI



Opal WE



Opal OC



Opal 56



Opal 57



Opal 58



Opal 59



Opal 60



Opal SL



Opal T



Opal AM-Y



Opal AM-V



Opal AM-R



OT



PT



BT



GT



T-GLASS



MY



Mlv



MP



Y



O



RED



G



W



CT-CL



CT-W



CT-A



CT-B



CT-R

Shade	W1	W2	W3	A1	A2	A3	A3.5	A4	rootA
Opaque Liner									
	W1O	W2O	W3O	A1O	A2O	A3O	A3.5O	A4O	rootAO
Cervical	—	—	—	—	—	A3B:2 AC:1	A3.5B:1 AC:1	A4B:1 AC:1	
	—	—	—	—	—	A3B:2 AC:1	A3.5B:1 AC:1	A4B:1 AC:1	AC
Margin		W1M:2 NM:1	W1M:1 NM:2		NM:1 A3M:1		A3M:2 rootAM:1	A3M:1 rootAM:2	
	W1M	W1M:2 NM:1	W1M:1 NM:2	NM	NM:1 A3M:1	A3M	A3M:2 rootAM:1	A3M:1 rootAM:2	rootAM
Opaque Dentin		OD-W1:2 OD-N:1	OD-W1:1 OD-N:2		OD-N:1 OD-A3:1		OD-A3:2 OD-rootA:1	OD-A3:1 OD-rootA:2	
	OD-W1	OD-W1:2 OD-N:1	OD-W1:1 OD-N:2	OD-N	OD-N:1 OD-A3:1	OD-A3	OD-A3:2 OD-rootA:1	OD-A3:1 OD-rootA:2	OD-rootA
Body									
	W1B	W2B	W3B	A1B	A2B	A3B	A3.5B	A4B	rootAB
Enamel		Opal 56:2 Opal 57:1	Opal 56:1 Opal 57:2				Opal 59:1 Opal 60:1		
	Opal 56	Opal 56:2 Opal 57:1	Opal 56:1 Opal 57:2	Opal 57	Opal 58	Opal 59	Opal 59:1 Opal 60:1	Opal 60	

Shade	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaque Liner											
	B1O	B2O	B3O	B4O	C1O	C2O	C3O	C4O	D2O	D3O	D4O
Cervical	—	—	B3B:1 BC:1		—	C2B:2 CC:1	C3B:1 CC:1				
	—	—	B3B:1 BC:1	BC	—	C2B:2 CC:1	C3B:1 CC:1	CC	DC		
Margin	NM:1 B2M:1		B2M:1 B4M:1		NM:1 C2M:1		C2M:1 C4M:1		NM:1 D3M:1		D3M:1 B4M:1
	NM:1 B2M:1	B2M	B2M:1 B4M:1	B4M	NM:1 C2M:1	C2M	C2M:1 C4M:1	C4M	NM:1 D3M:1	D3M	D3M:1 B4M:1
Opaque Dentin	OD-N:1 OD-B2:1		OD-B2:1 OD-B4:1		OD-N:1 OD-C2:1		OD-C2:1 OD-C4:1		OD-N:1 OD-D3:1		OD-D3:1 OD-B4:1
	OD-N:1 OD-B2:1	OD-B2	OD-B2:1 OD-B4:1	OD-B4	OD-N:1 OD-C2:1	OD-C2	OD-C2:1 OD-C4:1	OD-C4	OD-N:1 OD-D3:1	OD-D3	OD-D3:1 OD-B4:1
Body											
	B1B	B2B	B3B	B4B	C1B	C2B	C3B	C4B	D2B	D3B	D4B
Enamel											
	Opal 57	Opal 58	Opal 59	Opal 60	Opal 58	Opal 59	Opal 60	Opal 58	Opal 59		

Opal Effect								
	Opal T	Opal SL	Opal WE	Opal MI	Opal OC	Opal AM-R	Opal AM-Y	Opal AM-V

Enamel Effect						Cervical Trans					
	BT	OT	PT	GT	T-GLASS	Cervical Trans	CT-CL	CT-W	CT-A	CT-B	CT-R

Color Effect									Gum Color		
	MP	MY	Miv	RED	Y	O	G	W	Gum Color	Gum-LP	Gum-DP
						Correction					
	OM-Y	OM-LP	OM-DP	LPM	CLM	Correction			ADD-ON B	ADD-ON T	CPM-FINE

※印刷物のため実物とは色調が異なる場合があります。