

SPARE PARTS / ACCESSORIES

メルサーージュ エピック S

製品名	製品	製品番号		製品名	製品	製品番号	
フットコントローラー (エピックS用)		9966012		スケーラーハンドピース (ライト付)	 <small>[水] [熱]</small>	9966010	
水フィルターセット (2mブルー)		9966034		滅菌ケース	 <small>[熱]</small>	9966013	
水供給チューブ (2mブルー)		9966035		水供給用コネクタ		9966048	
ボトル400セット		9966037		ケース付水フィルター		9966049	
ボトル400用蓋 (エピックS用)		9966038		水フィルター		9966050	
ボトル400(蓋なし) (エピックS用)		9966039		チップ交換用レンチ (CR-10)	 <small>[水] [熱]</small>	9966052	
注水ポンプ (エピックS用)		9966040		Eチップ用レンチ	 <small>[水] [熱]</small>	9966053	
Oリング(細) (ボトル400用)		9966042		チップホルダー	 <small>[水] [熱]</small>	9966054	
Oリング(太) (ボトル400用)		9966043		チップカバー(L)	 <small>[水] [熱]</small>	9966055	
パッキン (エピックS 水ボトル蓋用)		9966045		チップカバー(S)	 <small>[水] [熱]</small>	9966056	
固定ネジ (S-V30/Vチップホルダー用)		9966065		電源コード		9966051	
チャック (S-V30/Vチップホルダー用)		9966066		チューブストッパー		9966046	
Oリング (S-E9/S-E12/S-V10用)		9966067		防塵カバー		9966044	
ナット (S-V10用)		9966068		スパナ (5×8) 2個入		9966057	
ナット (S-E9/S-E12用)		9966069		チューブ抑えリング (2個セット)		9966036	
スケーラー ハンドピースコード		9966011		ゴム足		9966047	
Oリング (スケーラーハンドピース コード用)		9966041					

[水] 熱水洗浄可能 [熱] 135℃までの温度でオートクレーブ滅菌可能